

Blood Bowl Templates - Instructions

This release contains scanned and homemade printable copies of the official templates used in the Blood Bowl miniatures game by GW.

Counter/KO/Dead/Reserves template one large A4 size and one small ½ A4 size

Scatter template

Sideline template

Throwing range ruler to be used with a field having about (the official) 29x29mm squares

All the templates are optimized for printing using an A4 color printer and then laminated. For all but the ruler this is pretty straightforward: 1) Print 2) Cut out 3) Laminate.

The ruler however comes in two parts since it is longer than a normal A4 paper and requires some assembling.

READ THE INSTRUCTIONS BEFORE STARTING!

1. Cut out the two parts of the ruler, but be sure to leave a bit of white before the LONG. This part will be used to glue the two parts together for stability when aligning for lamination.



2. Cut out the white hole in the checkered area on the QUICK-part. The hole is placed over your figure when measuring throws.



3. Glue the parts together using paper glue. You don't need much glue, just so the paper doesn't slide around.



4. a) If you have laminator plastic of A3 size you can cut a long stripe wide enough for the whole ruler to fit in, laminate and you are done.

A3 size plastic is however very expensive compared to A4 size plastic.

The following describes a way to use A4 size instead.

b) Cut an A4 sized plastic in almost half leaving one part slightly wider than the other, so you get two long stripes.



5. Put the QUICK-end in the widest stripe and the BOMB-part in the smaller stripe.



6. Now have the wide stripe overlap/enclose the smaller stripe. This will keep the plastic from curling inside your laminating machine.



7. Laminate, starting with the QUICK-end.

8. You are done; start playing blood bowl with your new templates.

